

SAT (450

TECHNICAL SPECIFICATIONS

System description	Benchtop, fully automatic, random access, open analyzer for clinical and immunoturbidimetric assays
ASSAY TYPES	- End point, Kinetic, Differential, Bichromatic, Fixed Time and Immunoturbidimetric - up to 999 tests, all user programmable
on-Board Reagents	Up to 72 reagents can be located into the refrigerated compartment and distribuited on: - 4 removable racks - 50/20/5 ml containers can be loaded - 8 extra positions for 5 ml containers
SAMPLE LOADING	 - 4 universal separate racks for continuous loading of samples, calibrators and controls thanks to a dedicated lid (new Main Cover Design) - each rack contains up to 15 cups and/or tubes 10-16 mm diam., 40-100 mm height - 8 extra positions for urgent samples, calibrators and controls - STAT sample execution is allowed at any time
THROUGHPUT	- Up to 440 tests/hour including I.S.E.
reading system	- Direct reading, dual channel photometer - 200 I minimum reading volume - 9 narrow band, automatically selected interferential filters 340, 380, 405, 492, 510, 546, 577, 620, 690 - halogen lamp 6V/10W - linear from 0.0001 to 4.200 Abs - high resolution (0.0005 Abs.)
Sampling process	One mechanical arm performs all sampling operations with: - level sensing system (capacitive) - new sampling probe with shock sensor - new reagent pre-warming system at 37°C - automatic probe washing - sample volume range: 2.0µl÷ 99µl(0.25 l incr.) - reagents volume range: 0µl÷ 350µl(1 l incr.) sample pre-dilution - sample post-dilution - sample post-concentration
REACTION PLATE	Holds 80 cuvettes divided into 4 racks with 20 cuvettes each and performs the following operations: - automatic cuvettes washing - continuous monitoring of cuvettes quality - incubation temperature control at 37°C ± 0.3
COMPUTER (minimum requirements)	CPU: i3 series 9 (or higher) MEMORY: Ram 8 Gb HARD DISK: 250 Gb MONITOR: 17" COMMUNICATION: USB and/or RS232 port HOST COMMUNICATION: RS232 or Ethernet SOFTWARE: Windows 10 or higher
OPTIONAL ACCESSORIES	- Direct potentiometry New I.S.E. module for ⁺ Na, ⁺ K , ⁻ Cl - positive barcode readers for reagents and samples (code 128, codebar, code 2 of 5 interleaved, code 39)
Installation requirements	Power Supply: 90 ÷ 250 VAC, 47 ÷ 63 Hz Electric Consumption: 386 Watt Water Consumption: 1.8 l/h
DIMENSIONS	- 107 x 68 x 53 cm (W x D x H)
WEIGHT	- 49 Kg
ENVIRONMENT	- Room Temperature: 18°C ÷ 30°C - relative Humidity: 20% ÷ 85%




AMS Srl

Instruments Factory, Sales & Customer Service Via E. Barsanti 17/A 00012 Guidonia - Rome (Italy)


Tel. +39 0774 354441 Fax +39 0774 578035 The SAT 450 is a fully automatic, random access analyzer, designed to support the modern laboratory in achieving improved cost efficiency, reduced rerun and manual handling. It provides quality results while following rigid quality system procedures and in strict compliance with regulatory agency requirements.

SIMPLICITY

& TECHNOLOGY



- real time presentation of samples and reagents status;
- dynamic monitoring of test execution;
- analyzer detailed status with error messages, incubation temperature control, START, STOP, and PAUSE functions.



- display of test reaction curves;
- display of test calibration curves;
- Levy-Jennings and Youden QC statistic.



- each instrument part can be singularly tested;
- remote diagnostic capability.

The **SAT 450** operating software runs under multitasking WINDOWS 10.

Its bright touchscreen software simplifies routine work and offers the operator a concise guide through the operating functions.



- all parameters are completely user programmable up to 999 methods can be memorized and all settings are user programmable;
- user selectable pre-dilution and re-dilution ratio, samples re-run option and reflex test;

AUTOMATIC DETECTION

- end point limit
- substrate depletion
- fixed time and kinetic reaction linearity

MORE FEATURES

- multilanguage files (user programmable)
- bi-directional communication compliant with ASTM and CLSI standard
- automatic transmission of patients data and results through serial port or Ethernet
- reagents data traceability and management
- Westgard rules management
- Sampling order with Priority management
- Possibility to compare current calibration with any previous calibration
- Preventive maintenance management
- Automatic startup and automatic shutdown procedures

CLINICAL CHEMISTRY

TOXICOLOGY

IMMUNOTURBIDIMETRY

PARAMETERS

Drugs of abuse

6-ACETYLMORPHINE

NORBUPRENORPHINE

BUPRENORPHINE

COTININE

AMPHETAMINE / METHAMPHETAMINE
BARBITURATES
BENZODIAZEPINES
CANNABINOIDS (THC)
COCAINE
ECSTASY (MDMA)
EDDP
METHADONE
OPIATES
ETHIL ALCOHOL
ETG (ETHYL GLUCORONIDE)
SYNTHETIC THC (SPICE K2)



Clinical chemistry

CC AMYLASE - Liquid Stable ALBUMIN - Liquid Stable BILIRUBIN Total - Liquid Stable BILIRUBIN Direct - Liquid Stable CALCIUM ARSENAZO - Liquid Stable CALCIUM OCP - Liquid Stable CHLORIDE - Liquid Stable TOTAL CHOLESTEROL - Liquid Stable HDL DIRECT - Liquid Stable LDL DIRECT - Liquid Stable HOMOCYSTEINE - Liquid Stable ENZYMATIC CREATININE - Liquid Stable CREATININE - Liquid Stable GAMMA-GT - Liquid Stable LDH - Liquid Stable PHOSPHOROUS UV - Liquid Stable GLUCOSE - Liquid Stable MAGNESIUM XB - Liquid Stable TOTAL PROTEIN - Liquid Stable PYROGALLOL PROTEIN - Liquid Stable TRIGLYCERIDES - Liquid Stable URIC ACID - Liquid Stable IRON FERENE - Liquid Stable GOT/AST - Liquid Stable GPT / ALT - Liquid Stable ALP - Liquid Stable UREA UV - Liquid Stable CHOLINESTERASE - Liquid Stable ENZYMATIC HbA1c - Liquid Stable GLICOTEST HbA1 - Liquid Stable CK-NAC - Liquid Stable ASO Latex - Liquid Stable CRP Latex – Liquid Stable RF Latex - Liquid Stable FERRITIN Latex - Liquid Stable

MICROALBUMIN Latex - Liquid Stable

IgA - Liquid Stable IgG - Liquid Stable

C3 - Liquid Stable

IgM - Liquid Stable

C4 - Liquid Stable

TRANSFERRIN - Liquid Stable

Method

Enzymatic CNPG3 Colorimetric BCG Colorimetric Jendrassik-Groff Colorimetric Jendrassik-Groff Colorimetric Colorimetric Colorimetric Enzymatic colorimetric Direct enzymatic Direct enzymatic Enzymatic Enzymatic colorimetric Jaffé modified Enzymatic Enzymatic DGKC Enzymatic Enzymatic colorimetric Colorimetric Colorimetric (Biurel) Colorimetric Enzymatic colorimetric Enzymatic colorimetric Colorimetric **Enzymatic IFCC-SCE** Enzymatic IFCC-SCE Enzymatic DGKC Enzymatic Enzymatic DGKC Direct Enzymatic Direct Latex Immunoturbidimetric Enzymatic IFCC **Immunoturbidimetric Immunoturbidimetric Immunoturbidimetric Immunoturbidimetric Immunoturbidimetric Immunoturbidimetric Immunoturbidimetric**

Immunoturbidimetric

Immunoturbidimetric

Immunoturbidimetric

COMPACT

SOLUTION

On board reagents

Up to 64 reagents on 4 separate racks additional 8 positions on board storage at 8°C ÷ 12°C barcode tracking of reagents (optional) possibility to close the system with dedicated reagents (optional)



Pratical design

Advanced technology and sophisticated system enginering allow easy access for maintenance purposes by simply lifting up the entire light-weight analytical chassis.



New I.S.E. Module (OPTIONAL)

direct measurement of Sodium, Potassium and Chloride easy access for reagents replacement maintenance-free electrodes Low Fluids Pack consumption





Simplicity, productivity, quality results

a single, heated level sensing arm with Shock Sensor minimizes mechanical complexity and rapidly moves to all positions for:

optimized sample and reagent mixing

up to 4 different reagents addition to each method sample pre-dilution

sample rerun with or without dilution

probe double washing between each sampling cycle inner and outer probe washing after every cycle



Cost Efficency

80 long life, 6 mm optical path,semi-disposable cuvettes (4 separate racks with 20 cuvettes each) selective cuvettes wash with 3 different liquids: acid, alkaline and neutral low water consumption: 1.8 l/h 200 l minimum reaction volume cuvettes quality automatic check



Flexibility

4 universal racks for cups and primary tubes continuous loading of up to 15 samples on each rack (serum, plasma, urine) thanks to dedicated lid additional 8 positions

STAT samples loading allowed at any time standards and controls can be placed on any avaible rack with no limitations positive identification of samples (optional)