medline

scientific limited

HEATING / STIRRER MANTLES

OPERATING MANUAL

HEATING MANTLE INSTRUCTIONS

On opening the packaging please check that the Mantle is undamaged and that it is the correct voltage.

Construction:

The mantles are constructed from a coated steel body with a soft glass fibre mantle. Insulation is manufactured from compressed silica and the "ON" status is shown by a power indicator lamp. In the EAM series there is an additional stainless steel cover for both three and six place models.

The DM and CM series also include a steel temperature probe

Operation:

Place the pre filled flask into the mantle ensuring the flask is the correct size to give good thermal conductivity. Smaller flasks can be placed in to mantles if supported but will be limited to mainly bottom surface heating. Plug the unit into the power source or switch on the power source if already connected.

In the case of the C series, the mantle will now heat to approximately 450c and

regulation can be provided by an external regulator.

In the case of the E series, temperature setting is regulated by the analogue dial and the required level is determined by using a thermometer (not supplied). This also applies to the EAM series with individual settings for each position.

For the CM series the temperature is reached by setting the dial to the required level and checked with a thermometer if required. Monitoring and stability are achieved with the K type probe which is plugged into the rear of the unit.

For the DM series the temp level is set by use of the digital potentiometer and monitored by the K type probe.

The ES series Stirrer Mantles have an additional speed (RPM) setting.

During the heating stage on the CM, E, EAM, ES series the yellow "Heating" light will come on periodically. The frequency of the heating mode will depend on the temperature setting.

Turn off the heating controllers before removing the flasks and when the process is

finished the unit must be turned off and the plug removed from the power source.

After using the instrument for the first time the white heating element will become discolored. This is normal due to the high heating temperatures and will not affect the performance of the product.

Important Safety Precautions:

Do not use the mantle with metal or conductive flasks.

Do not place or remove the mantle when switched on.

Do not touch the element or flask when in use.

Avoid spilling any liquid on the flask or mantle.

In Case of Accident

Should spillage occur onto the heating mantle, immediately turn off the power at the mains and remove the IEC type power lead from the mantle body.

Remove any broken glass from the mantle element when cool using protective hand covering to avoid dangers of acids, contamination and cuts.

Place the mantle in an oven set to below 60C° or a dry, warm place until dry.

Before reuse check the functioning of the mantle with a flask half filled with water. Should the instrument malfunction in anyway please contact the supplier for details of servicing the unit or replacement.

The Medline Range Consists Of:

C – Series:	Uncontrolled
CM – Series:	Controlled (with temperature controller)
DM – Series:	Digital (with temperature controller)
E – Series:	Controlled
ES – Series:	Stirrer
EAM – Series:	Multi Position (3 or 6)

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Medline Scientific Limited